



## Advanced Nanoscale Characterization with diInnova™ *Phase Imaging of Polymer Materials*

By: Thomas Müller, PhD, Applications Scientist  
Veeco Instruments Inc.

### INTRODUCTION

TappingMode™ imaging has proved to be the most versatile mode of atomic force microscopy (AFM) in ambient conditions where the presence of a fluid layer (condensed water vapor and other contaminants) severely limits the applicability of both, contact mode and non-contact techniques. Overcoming the challenges posed by friction, adhesion, and other issues, TappingMode has provided a means of greatly extending AFM applications, particularly regarding high resolution imaging of soft samples.

Phase Imaging is an important extension of TappingMode imaging. By mapping out the phase of the oscillating cantilever, phase imaging goes beyond simple topographical mapping. Being sensitive to variations

in adhesion and viscoelasticity, phase imaging can provide information about sample composition and microphase separation. A key strength of phase imaging is its ability to provide such information at the high resolution achievable with TappingMode.

First class high resolution performance is critical for taking full advantage of TappingMode phase imaging. With its stable, low-drift platform, ultra-low noise closed loop scan control, and excellent force control the new diInnova SPM (see Figure 1) is an ideal instrument for high resolution imaging of delicate samples. In addition, the Innova combines this outstanding core performance with generous data acquisition bandwidth and facile signal access, thus enabling a wide range of demanding research applications.



Figure 1: The new Innova

## APPLICATION EXAMPLES

Phase imaging can reveal microphase separations occurring in block copolymers. Obtaining this information with alternative techniques involves complications, such as chemical staining for TEM. With TappingMode phase imaging, AFM can provide the visualization of the microphase separation pattern directly from images obtained in ambient conditions of an untreated thin film. Figure 2 shows topography and phase images of a PS-*b*-PB-*b*-PS triblock copolymer (PS, polystyrene; PB, polybutadiene)<sup>1</sup>. Both channels clearly show the expected worm-like microphase separation pattern. The microphase domains exhibit a width of ~ 35nm. This can be seen most clearly in the 2-D Fourier Transform shown in Figure 3. Note that the intensity maximum is perfectly circular, as it should be given the overall isotropy with no preferred block orientation and no dependence of block width on azimuthal angle. These images were acquired using closed-loop scan control, ensuring calibrated, undistorted measurements.

For a phase image to reflect differences in viscoelasticity or modulus of a

material, the AFM probe needs to penetrate the material. More precisely, the probe needs to penetrate sufficiently far such that the tip-sample interactions are influenced by material properties from the layer of interest. In the case of a PS-*b*-PB-*b*-PS film, an amorphous, PB-enriched top-layer is usually present. Thus, the combination of a soft cantilever (e.g., FESP,  $k \sim 2\text{--}5\text{N/m}$ ) with very light tapping conditions would fail to uncover the microphase separation pattern. Images such as those shown in Figure 2 are usually obtained with fairly hard tapping conditions, that is, fairly high ratios of free amplitude to amplitude setpoint. On Innova, probe tuning at moderate amplitude (input gain setting ~ 8 or 10) and significant increase of the cantilever drive amplitude upon engaging will yield the desired result.

As phase imaging is sensitive to local variations in mechanical properties, it can afford an efficient means for mapping out the distribution of components in composite samples. Figure 4 shows topography and phase images of a cross-sectioned multilayer polyethylene sample, composed of alternate high and low density layers.

The topography image is dominated by the large scale, low frequency height undulation that has apparently resulted from cross-sectioning by cryo-microtoming. The phase image has a different appearance, clearly providing complementary information. It is dominated by an alternating set of stripes, obviously representing the sought after alternation in material properties and thus component layers. In addition, the phase image reveals topographic fine features that are much less apparent in the height image. In particular, small droplets can be discerned with well-defined phase contrast indicating distinct local mechanical properties. The formation and coalescence of small droplets on microtomed polyethylene samples indicate an ageing surface. Note that the droplets are not distributed randomly. Rather, some of them appear to form along lines, presumably small scratches imparted on the sample by the microtoming process.

In both, Figures 2 and 4, the phase image brings out the material contrast most clearly and separates it from information about large scale topographic information. However,




Figure 2: Topography (left) and phase image (right) of PS-*b*-PB-*b*-PS triblock copolymer. Sufficiently hard tapping conditions have ensured probe penetration into the subsurface layer, where a wormlike microphase separation pattern is present as can be seen clearly in both channels. Image size 2.0 $\mu\text{m}$ . Closed loop active.

<sup>1</sup>Samples courtesy of N. Erina

in both cases, the material contrast seems to be also partly contained in the topographic image. The parts of the sample appearing with bright contrast in the phase image appear raised in the topographic image. This can be rationalized as a reflection of

material stiffness. Under hard tapping conditions, the probe penetrates into the material. The regions appearing bright in the phase image are stiffer, leading to less probe penetration and thus a raised topography relative to the softer parts at high force. A

more complete explanation of this effect has to include the nature of the feedback in tapping mode. Strong positive phase shifts in hard tapping conditions indicate significant up-shifts in the resonance frequency of the cantilever – while drive frequency and amplitude setpoint chosen for feedback remain constant and were chosen at the resonance frequency of the free cantilever. In the regions with the most positive phase shift, the resonance frequency effectively shifts furthest away from the drive frequency. Being furthest off-resonance, the cantilever is being driven less efficiently (while the amplitude setpoint remains unchanged), leading effectively to lighter tapping, which contributes to the appearance of a raised profile.

Aside from compositional mapping and the visualization of microphase separations, phase imaging can aid in the detection of fine structures. In the case of the MLPE sample shown in Figure 4, interesting fine structure can be observed in higher resolution images of the layer boundaries when slightly lighter tapping conditions are employed. As can be seen in Figure 5, hair-like structures appear to extend from the boundary into the layer appearing with darker phase




Figure 3: 2-D Fourier transform of the phase data shown in Figure 2. The ring-shaped intensity maximum indicates that the phase separation pattern is isotropic with a well-defined repeat distance of  $r=35\text{nm}$  as indicated at the bottom of the dialog and as expected for this triblock copolymer.




Figure 4: Topography (left) and phase image (right) of a cryo-microtomed multilayer polyethylene sample. While topography is dominated by large-scale undulations, phase provides a clean view of the layered structure. Additional fine structure shows the presence of small droplets. Image size  $35\mu\text{m}$ . Closed loop active.




Figure 5: Phase image of a cryo-microtomed multilayer polyethylene sample. Hair-like fine structure can be seen near the layer-interface. Image size 5µm. Closed loop active.

contrast. Closer inspection of Figure 5 reveals lamellar structures throughout the lower density (darker phase) component visible in the right half of the image and a gradual loss of alignment with increasing distance from the interface.

Figure 6 shows AFM images of an oriented film of isotactic polypropylene, also known as the microporous membrane Celgard. Both, topography and phase clearly show the pattern of oriented fibrillar structures that is characteristic of this sample. With

the overall height scale (~ 200nm) dominated by large variations, finer structures are not evident in the topography data. In contrast, the phase image exhibits very clear and well-defined additional features. Fine lamellar structures (only ~ 20nm wide) are seen to be present in between the rows of fibrils. The lamellar structures are seen to be oriented perpendicular to the larger fibrillar structures. As the phase signal is sensitive to deviations of the oscillation amplitude from the amplitude setpoint, it can serve as an edge detection technique and thus highlights such fine structures that are easily overlooked in the topography channel. The accurate imaging of fine structures on this delicate sample benefits greatly from Innova's combination of excellent force control and low noise closed-loop scan control.

The appearance of fine structure in phase images complements the sensitivity to material properties. By identifying components in composite samples, the appearance of fine structure in phase images aids in compositional imaging, see (AN103) *Phase Imaging of Polymer Materials with diCaliber*. Even higher resolution phase images can reveal length scales

associated with the self-assembly of individual molecules in monolayer films and their relationship to the substrate. Figure 7 shows topography and phase images of a self-assembled monolayer of  $C_{60}H_{122}$  alkane molecules on a substrate of highly-oriented pyrolytic graphite (HOPG).

Both images clearly show the presence of patches (domains) separated by sharp (mostly straight) borders. Closer inspection reveals a pattern of straight lines within each domain. This lamellar pattern is much more obvious in the higher resolution phase image shown in Figure 8. The lamellae are seen to have a well-defined spacing as well as preferred directions. This is confirmed in Figure 9, which shows the Fourier Transform of the phase image shown in Figure 7. The Fourier Transform shows very clearly the hexagonal symmetry of the self-assembly pattern that reflects the hexagonal symmetry of the underlying graphite substrate. When forming self-assembled monolayers,  $C_{60}H_{122}$  alkane retains a fixed relationship to the high symmetry axes of graphite.

The special relationship of the




Figure 6: Topography (left) and phase image (right) of Celgard. While oriented fibrillar structures are evident in topography, the phase image additionally reveals lamellar fine structure. Image size 2.5µm. Closed loop active.




Figure 7: Topography (left) and phase image (right) of a  $C_{60}H_{122}$  monolayer self-assembled on graphite. Both images clearly show self-assembled domains, each of which is composed of parallel lines. Image size  $1.5\mu m$ . Closed loop active.




Figure 8: Phase image of a  $C_{60}H_{122}$  monolayer self-assembled on graphite clearly showing lamellar fine structure associated with the self-assembly. Image size  $390nm$ . Closed loop active.

adsorbate structure with the graphite substrate is consistent with the assumption, that the layer probed here is actually the molecular layer that is in direct contact with the substrate. The question arises because the preparation of the  $C_{60}H_{122}$  alkane sample cannot be assumed to result in a single molecular layer. Indeed, judicious choice and excellent control of tapping forces is required to reveal the structures shown here. The probe has to penetrate through partially disordered, soft adsorbate multilayers without destroying that “first”




Figure 9: 2-D Fourier Transform of the phase image shown in Figure 7. The hexagonal symmetry is clearly visible. As shown near the bottom, the periodicity is measured to be  $r=7.5nm$ .

molecular layer that is in direct contact with the substrate and is thereby subject to additional stabilization.

As indicated in Figure 9, the spatial periodicity is seen to be about  $7.5nm$ .  $C_{60}H_{122}$  alkane is known to self-assemble on graphite such that each molecule assumes an extended all-trans conformation with its backbone parallel to the substrate and perpendicular to the lamella axis. Therefore, the lamella width equals the length of a single  $C_{60}H_{122}$  alkane molecule, which is about  $7.5nm$ . For

more on probing the self-assembly of alkanes and their derivatives on graphite see (AN85) *Scanning Tunneling Microscopy: A Tool for Studying Self-Assembly and Model Systems for Molecular Devices*.

The interrogation of self-assembly in monolayers on flat substrates depends critically on good force control for nondestructive imaging, in conjunction with high stability and low stage drift to enable the required high resolution performance. Obviously, the analysis presented above also requires correct scanner calibration. Closed-loop scan control can ensure correct calibration but is often associated with excessive noise levels. Not so on Innova. In fact, all images presented in this application note (including the  $390nm$  image shown in Figure 8) were acquired on Innova with closed-loop control active, thus ensuring accurate measurements.

#### SUMMARY

With TappingMode phase imaging, the Innova system can efficiently and nondestructively map variations in sample properties at the highest resolution. As TappingMode is often the preferred imaging mode for delicate samples, phase imaging can complement other modes such as

force modulation and lateral force microscopy, often with superior image detail. Phase imaging applications include the characterization of composite materials, mapping of variations in adhesion and viscoelasticity, and identification of surface contamination. Highest resolution phase images open the door to studies of molecular self-assembly. The combination of excellent high resolution performance and TappingMode phase imaging makes Innova a powerful tool for the study of material properties at the nanometer scale.

## ABOUT THE IMAGING TECHNIQUES

### Conventional Approaches

Two conventional AFM scanning modes – contact mode and noncontact mode have been used for some time with varying success. Each has its limitations, particularly for imaging delicate samples in ambient conditions, see also (AN04) *TappingMode Imaging Applications and Technology*.

Contact mode AFM represents the simplest imaging technique. The sample is simply moved laterally relative to the probe such that the probe is dragged across the surface. While this technique has been successful for many samples, it is subject to serious drawbacks. In essence, the dragging motion of the probe combined with adhesive tip-surface forces can lead to substantial damage to probe and sample, creating artifacts in the image and often degrading the resolution severely.

Under ambient air conditions, most surfaces are covered by a fluid layer, composed of water and other contaminants, which is typically several nanometers thick. An AFM tip touching this layer will cause a meniscus to form and surface tension will pull the tip onto the surface. Additional adhesive forces can arise from trapped electrostatic charges (see Figure 10). Tip-sample forces are thus larger than the cantilever deflection would seem to indicate. Correspondingly, lateral motion during contact mode imaging is

associated with larger lateral forces, resulting in severe tip or sample damage or involuntary displacement of weakly bound surface adsorbates. Capillary forces can be eliminated by completely submersing the sample and probe tip in liquid. However, sample surfaces are often either less robust in liquid (e.g., adsorbates are more weakly bound) or are not compatible with a liquid environment at all.

Non-contact mode represents an attempt to overcome the deleterious tip-surface forces associated with contact mode. Unfortunately, ambient air conditions are rarely conducive to non-contact AFM imaging.

Non-contact mode is based upon the detection of weak attractive van der Waals forces that exist between tip and sample a few nanometers above the surface. It is important to note that the fluid layer present in ambient conditions partially shields these forces and occupies a large fraction of their useful range, i.e., when compared to operation in ultra-high vacuum, where no fluid layer is present. Due to the limited range of any residual van der Waals forces, their detection necessitates very small




Figure 10: In contact mode AFM, electrostatic and/or surface tension forces from the adsorbed fluid layer lead to destructive lateral shear forces. From application note AN04, *TappingMode Imaging Applications and Technology*.

oscillation amplitudes. At the same time, a cantilever operated at very small amplitudes is easily trapped inside the fluid layer, once it touches. In the ideal case, operation remains outside the fluid layer – and therefore several nanometers away from the sample surface. The consequences are substantially degraded resolution (as compared with TappingMode) and an inability to map out variations in

local mechanical properties (see the preceding sections of this application note for many examples of probing local mechanical properties with TappingMode phase imaging). In practice, the probe is frequently drawn onto the sample surface and remains trapped there while the scanning motion proceeds, leading to unusable data and sample damage similar to that caused by contact mode.

### TappingMode Imaging in Air

Veeco's patented TappingMode imaging overcomes the limitations of the conventional scanning modes by alternately placing the tip in contact with the surface to provide high resolution and then lifting the tip off the surface to avoid dragging the tip across the sample. TappingMode imaging is implemented in ambient air by oscillating the cantilever at or near its fundamental flexural resonance, usually a frequency in the range of 50 to 500 kHz. "Free air" amplitudes are typically greater than 20nm. During the engage process, the tip-sample separation is reduced until the tip begins to interact with the surface. The interaction with the surface ("tapping") leads to energy loss and a reduced oscillation amplitude, see Figure 11 and (AN04) *TappingMode Imaging Applications and Technology*. Deviations of the amplitude from a setpoint value serve as error signal in the feedback loop that drives the Z-motion to track surface features. In phase imaging, the phase lag of the cantilever oscillation, relative to the drive signal, is simultaneously monitored with topography data. As the phase lag is influenced by energy dissipation experienced during the




Figure 11: TappingMode cantilever oscillation amplitude in free air and during scanning. From application note AN04, *TappingMode Imaging Applications and Technology*.




Figure 12: Selecting drive frequency and phase in preparation for TappingMode imaging using the autotune function in the Cantilever Tune dialog.

oscillation cycle, it is very sensitive to material properties such as adhesion, viscoelasticity, and modulus. Phase imaging can also act as a real-time contrast enhancement technique. Because phase imaging highlights edges and is not affected by large-scale height differences, it provides for clearer observation of fine features, such as grain edges, which can be obscured by rough topography.

TappingMode prevents the tip from sticking to the surface and causing damage during scanning. In contrast to non-contact mode, the oscillation amplitude and energy is sufficient to overcome tip-sample adhesion in every oscillation cycle. In addition, TappingMode provides a large linear operating range, i.e., linear dependence of cantilever amplitude on tip-sample separation, making the feedback system highly stable and allowing routine reproducible sample measurements at high resolution.

Selection of the optimal cantilever oscillation frequency is a critical step in preparing for TappingMode imaging. On Innova this is accomplished in the

Cantilever Tuning dialog (see Figure 12). Choice of drive frequency and adjustment of lock-in parameters for proper phase imaging are all taken care of by the AutoTune function. The user chooses the desired cantilever oscillation amplitude by selecting a target amplitude and lock-in sensitivity ("input gain") factor. As shown in Figure 12, both, amplitude and phase are displayed in the Cantilever Tuning dialog.

## REFERENCES

Q. Zhong, D. Inniss, K. Kjoller, and V.B. Elings, *Surface Science* 290, L688-92 (1993)

V. V. Tsukruk, *Rubber Chemistry and Technology*, 70(3), 430-67 (1997)

S. S. Sheiko, "Imaging of Polymers Using Scanning Force Microscopy: From Superstructures to Individual Molecules," *Advances in Polymer Science*, v. 151, p. 61-174, Springer, Berlin 2000

S. N. Magonov, "AFM in Analysis of Polymers," *Encyclopedia of Analytical Chemistry*, R. A. Meyers, Ed., p. 7432-91, John Wiley & Sons Ltd, Chichester, 2000.

D. A. Chernoff, S. N. Magonov, "Atomic Force Microscopy", Chapter 19 in "Comprehensive Desk Reference of Polymer Characterization and Analysis", R. Brady, ed., Oxford University Press, New York. 2003.



### WORLDWIDE CUSTOMER SUPPORT FROM THE INDUSTRY LEADER

Bruker Corporation is a leading provider of high-performance scientific instruments and solutions for molecular and materials research, as well as for industrial and applied analysis. For more information, visit [www.bruker.com](http://www.bruker.com), email [productinfo@bruker-nano.com](mailto:productinfo@bruker-nano.com), or call +1.805.967.1400/800.873.9750.