

## Application Note # AFM-02

# Investigation of Solar Cells

Solar cells, or photovoltaic cells, are used to convert sunlight into electrical power. As traditional power sources grow scarce, other forms of producing electrical power are gaining firm footing in the power supply mix. Solar cells are already widely used in a variety of applications - from spacecraft, to small portable devices, to farm installations, to roadway signs. As energy prices increase, public demand for solar power has surged. In order to meet the longevity, yield, and price requirements of consumers and industry, public and private sector research has increased dramatically.

While traditional tools are helpful to investigate and improve solar cells, AFM/SPM offers metrology, topography & roughness analysis at much higher resolution than with optical techniques. In addition AFM/SPM techniques can be used to measure electronic properties of samples, such as work function.




Fig. 1: Close up of the N8 TITANOS sample stage, showing the confocal microscope, a video microscope and NANOS (AFM/SPM scanning head). The sample stage accommodates samples up to 300 mm x 300 mm. Every position on the sample is available for rapid inspection.

## Measurements on solar cell components




Fig. 2:  
 a) Confocal microscope image of a thin film solar cell, image size:  $158\ \mu\text{m} \times 158\ \mu\text{m}$ , Z-height:  $6\ \mu\text{m}$   
 b) Intermittent contact mode AFM scan of ZnO, image size:  $10\ \mu\text{m} \times 10\ \mu\text{m}$ , Z-height:  $673\ \text{nm}$   
 c) Contact mode AFM scan of  $\text{SnO}_2$ , image size:  $3\ \mu\text{m} \times 3\ \mu\text{m}$ , Z-height:  $229\ \text{nm}$

## Technical details

All measurements were performed on an N8 TITANOS built for the FZ Jülich, Institute for Energy Research (Photovoltaics), and modified slightly to include a high-performance confocal microscope instead of the standard optical microscope. The N8 TITANOS is a high-resolution, large sample AFM that can measure with atomic-scale resolution anywhere on a 300 mm wafer, with precise repositioning under two or more measurement stations (AFM, confocal microscope etc).

The thin film solar cells were kindly supplied by the FZ Jülich, IER. Research into the Transparent Conductive Oxide (TCO) is fundamental for improving performance of Silicon thin film solar cells. Si thin film solar cells promise higher flexibility than present wafer-based cells, whilst significantly reducing manufacturing costs.

## Samples

All samples courtesy of the Institute of Energy Research (Photovoltaics), FZ Jülich.

## Authors

B. L. Holmes, A. Thöne  
 Bruker Nano, Herzogenrath, Germany

## Laser lines for interconnecting solar cell modules




Fig. 3:  
 a) AFM scan, image size:  $30\ \mu\text{m} \times 30\ \mu\text{m}$ , Z-height:  $1973\ \text{nm}$   
 b) Confocal image, image size:  $158\ \mu\text{m} \times 158\ \mu\text{m}$ , Z-height:  $8.194\ \mu\text{m}$

All configurations and specifications are subject to change without notice.  
 © 2010 Bruker Nano. Printed in Germany